

# P501/P601/P602/P603/P604: PEPTIDE MAPPING & PROTEIN THERAPEUTICS CHARACTERIZATION

## APPLICATION

- Comprehensive characterization of purified proteins to be developed for therapeutics
- Protein sequence confirmation and PTM mapping
- Comparative analysis of biosimilars with reference standards
- Evaluation of purified or partially purified protein samples including therapeutics products at the early stage of development
- Confirmation of full amino acid sequence and PTM mapping;
- N-linked glycosylation sites mapping and profiling of N-linked glycans of pure mAbs

## SUITABLE SAMPLE TYPE

- Purified protein, in solution or dried

## MINIMUM SAMPLE REQUIREMENT

- Purity >95%
- $\geq 10$   $\mu\text{g}$ /sample is recommended for optimal results
- Detection limit:  $\geq 0.1$   $\mu\text{g}$
- *Please contact us for any specific samples or requirements*

## SAMPLE PREPARATION

- Store the samples at  $-20^{\circ}\text{C}$  or  $-80^{\circ}\text{C}$  prior to shipping
- Package the samples with dry-ice or chemical blue-ice bag when shipping

## GENERAL SAMPLE GUIDELINES FOR MASS SPECTROMETRY ANALYSIS

- ✓ **Protein Identification:** Any biological sample; detection limit  $\geq 1$  ng/purified protein,  $\geq 0.1$   $\mu\text{g}$ /mixture of proteins
- ✓ **Protein Profiling (Quantitative Proteomics):** Any biological sample; (e.g., tissue, cells, protein lysates, serum/plasma, etc.) with minimum requirements of  $>100$   $\mu\text{g}$ /sample of total protein;  $>5$  mg/sample for tissue; or  $>1$  million cells/sample for cell pellets
- ✓ **PTM Analysis:** Phosphorylation, Acetylation, Methylation, Ubiquitination, etc. with minimum sample requirement of target protein  $\geq 0.1$   $\mu\text{g}$ /sample
- ✓ **Intact Protein MW Analysis:**  $\geq 5$   $\mu\text{g}$ /sample of any purified protein  $\leq 180$  kDa
- ✓ **Peptide mapping:**  $\geq 20$   $\mu\text{g}$ /sample of any purified protein with a purity  $>90\%$  for 100% sequence coverage